## Manual • ENGLISH

oxyhero push is a nasal dilator which softly widens the nasal wings allowing you to breathe more effectively through the nose during sleep and also while doing sports.

During physical activity nasal respiration will enhance breathing performance (lung capacity) by up to 30% - especially with regards to popular endurance sports such as running and cycling.

During sleep, unimpeded nasal respiration ensures that less turbulence is created in the upper airways, therefore reducing the amount and intensity of vibrations affecting pharyngeal tissue. As a result

thereof, snoring will become less of an issue. Moreover, you will also improve your body's regeneration capacity due to a rise in blood oxygen saturation by up to 15%.

There are also mini pollen filters available for the oxyhero push, which can be pushed into the nasal dilator's cylinders. They are made of refined wool (sheep wool) that provides them with a natural self-cleaning function – dirt particles, particulate matter or pollen are filtered and the risk of getting an allergic reaction is therefore minimised. The filters can be ordered separately at: [www.oxyhero.com](http://www.oxyhero.com)

### This is how you best use the oxyhero push:

**1** First clean your nasal passages. Nasal douches can be very helpful for that purpose. Consider using a nasal irrigation device such as Nose Buddy.

**2** When making use of the nasal dilator for the first time, standing in front of a mirror will help you get accustomed to it much quicker.

**3** Remove the oxyhero push size M from the box. The nasal dilator consists of two cylinders, which are conjoined via a connecting bridge. On the bottom of each cylinder there are two soft, thin blades. These blades are used for better grip in the nose, because they create a „plug“ effect. If you do not need these blades or find them unpleasant, you can simply shorten or cut them off with scissors.

**4** When inserting, hold the product so that the cylindrical tubes are facing up (the horseshoe-shaped connecting piece is facing down). Now clasp the blades on each tube with your fingers so that they are pressed against the cylinder and the cylindrical tubes face upwards. Now slide the oxyhero push with slight pressure into your nose. The nasal

dilator should sit firmly and comfortably - the blades give it extra support as they smoothly spread the nose outwards (like a plug).

**5** If the device does not sit tight enough, do step 4 again using oxyhero push in size L. If it sits too tight, please consider size S instead. You can also moisten the nasal dilator with water prior to use, which will make it easier for you to insert the device. If you suffer from skin irritation, you can further apply nasal balm.

**6** If you have a very small and narrow nose, it may be that you will not be able to insert the oxyhero push as deeply as it is usually intended (in this case consider cutting off the blades). The nasal dilator should not cause any pain while wearing it. Needless to say, you will have a slight sensation of a foreign body at the beginning, but which will disappear after a short period of getting accustomed to it.

**7** To remove the nasal dilator, simply pull the oxyhero push gently out of the nose.

### Cleaning instruction:

The oxyhero push should be thoroughly cleaned with warm water after each application. Afterwards use a clean towel to dry it. Refrain from strong and alkaline-based detergents while cleaning the device, as they affect the material and should not get in contact with the sensitive nasal mucosae. The oxyhero push should be stored in a clean and dry place.

### Caution:

small part - please keep safe from infants (risk of suffocation)! Not suitable for children under 12 years - use only under parental supervision!

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