
Instructions for Use and Care

Indications for Use

The ZQuiet® anti-snoring device is intended for the treatment of nighttime snoring in adults.



Why ZQuiet® Works

Snoring is a vibration of soft tissue against the back of the throat. As we sleep, the tongue relaxes and partially closes the airway causing vibration, or snoring, as air passes in and out. The ZQuiet® device repositions your lower jaw slightly forward, widening the airway to improve airflow and eliminate the snoring noise.


Instructions for Use

- Take a few moments to familiarize yourself with the various features of ZQuiet®, (see Diagram 1) Clean each device prior to initial use (follow instructions below).
- Identify Device #1 and Device #2 by the number on the underside of the upper tray, (see Diagram 1)

- Start with Device #1 your first night. Before sleep, place ZQuiet® in your mouth.
- When inserting ZQuiet®, hold the mouthpiece in the folded position with the upper tray facing up. Insert ZQuiet® with the hinge towards the rear of your mouth (see Diagram 2) allowing your teeth to settle into the upper and lower trays. You should notice that your lower jaw moves slightly forward of its normal position.
- If your snoring noise stops, continue using Device #1. If, after three nights using Device #1, your snoring noise has not stopped, use Device #2, which advances your lower jaw slightly more than Device #1.

- Upon waking, gently remove your ZQuiet® device and follow the recommended cleaning and storage instructions.

Comfort Tips

You may experience an adjustment period with ZQuiet®. Your jaw muscles may experience temporary soreness, like the soreness that may happen to muscles when trying a new exercise.

- First-time users can soak the device in hot tap water for one minute to slightly soften the device prior to inserting in mouth.
- Initially, try wearing the device for short periods of time during waking hours. Insert ZQuiet® for 20-30 minutes while reading or watching TV.
- If your jaw muscles are particularly sore, take a night off.
- You can trim ZQuiet® with small scissors to improve fit and comfort. After trimming, use a fingernail file or emery board to smooth any rough edges.

Cleaning Instructions

Clean your ZQuiet®, each morning to minimize harmful germs and bacteria as follows:

- Using a toothbrush, scrub all surfaces with very warm water and a small amount of mild soap or detergent.
- Rinse thoroughly with clean warm water to remove all soap.
- After rinsing, shake off all excess water and pat dry with a soft clean towel.
- Store ZQuiet®, at room temperature in the supplied case.

Important: Do not soak ZQuiet® in mouthwash, especially if it contains alcohol.

Possible Side Effects

You may experience some of the following reactions to ZQuiet®. These are normal and should disappear with continued use:

- You may feel pressure against your front teeth.
- You may experience a gagging sensation or an urge to cough.
- You may experience increased salivation.
- ZQuiet® may be inadvertently removed from the mouth during sleep.
- Upon waking in the morning, the lower jaw may continue to feel more forward than usual or you may experience a sensation of bite change. This should go away shortly after the next hour.

If you experience sharp or persistent pain lasting more than two hours after use, you should discontinue use and contact your dentist.

Contraindications

- Do not use ZQuiet® if you have been diagnosed with Central Sleep Apnea.
- Do not use ZQuiet® if you have been diagnosed with severe respiratory disorders.
- Do not use ZQuiet® if you have loose teeth, loose dental work, orthodontia, full dentures or advanced periodontal disease.
- Do not use ZQuiet® if you have been diagnosed with temporomandibular joint disorders (TMD), experience clicking of the jaw and/or have pain in the jaw joints.
- Do not use ZQuiet® if you are under the age of 18 years.

Warnings

Usage of ZQuiet® may cause:

- Tooth movement or changes in dental occlusion.
- Gingival or dental soreness.
- Pain or soreness to the temporomandibular joint.
- Obstruction of oral breathing.
- Excessive salivation.

Manufactured by:

ZQuiet®
5247 Shelburne Rd., Suite 204
Shelburne, VT 05482
U.S.A.

Copyright © 2014 Sleeping Well, LLC. All Rights Reserved
ZQuiet® is a Registered Trademark of Sleeping Well, LLC.
U.S. Patent #8166976, Foreign Patents Pending
ZQFA-IFU 04.14

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Copyright © 2014 Sleeping Well, LLC. Alle Rechte vorbehalten
Copyright © 2014 Sleeping Well, LLC. Tous droits réservés.
Copyright © 2014 Sleeping Well, LLC. Alle rechten voorbehouden
Copyright © 2014 Sleeping Well, LLC. 保留一切權利