somnipax
guard standard

&

somnipax
guard large

**Gebrauchsinformation/Manual/Notice d'utilisation/Indicazioni
d'uso/Instrucciones para el uso/Handleiding/Käyttöohje**

**Lieferumfang/Delivery includes/Contenu de la livraison/Ambito Fornitura /
Alcance de la entrega/Omvang van de levering/Pakkauksen sisältö**

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Manual • ENGLISH

somnipax guard (standard / large)

1. What is the somnipax guard (standard / large) mandibular advancement device and what is it used for?

somnipax guard (standard / large) is a two-piece jaw splint made of thermoplastic material. It's worn inside your mouth to prevent simple snoring (also known as primary snoring) and mild to moderate obstructive sleep apnoea.

The mouth guard is available in 2 different versions: somnipax guard standard is suitable for small to medium-sized jaws, and somnipax guard large is designed for large jaws. Please use our instructions for working out the appropriate size before starting the adjustment process.

somnipax guard (standard / large) is a class I medical device.

2. Intended use and scope of application

Intended use

This mandibular advancement device is used to prevent snoring and interruptions in breathing during the night. Snoring occurs when tissue in the upper airways slackens and the airway is narrowed. The air flows more rapidly at the points of constriction and causes the surrounding tissue to vibrate. This creates the sound we recognise as snoring. In the case of nighttime interruptions in breathing, the slackened tissue not only narrows the airway, but temporarily blocks it, either completely or partially. Breathing comes to a brief halt.

Scope of application

The mandibular advancement device is used to prevent the noise produced by simple positional snoring (also known as tongue snoring) and to prevent interruptions in breathing caused by positional, mild, or moderate obstructive sleep apnoea (also: unhealthy snoring). As a general rule, the mandibular advancement device is used alone. However, it can also be used in conjunction with positional therapy products (such as an anti-snoring T-shirt or an electronic positioning belt). somnipax guard (standard / large) must only be used for these purposes.

The mandibular advancement device should not be used if any of the following apply:

- If the mouthpiece cannot be firmly positioned against the upper and lower jaw in a way that is sufficiently stable. E.g. if a full or partial removable denture is worn; or if there are fewer than 10 teeth seated firmly per jaw (dental fillings do not present any impediment).
- If braces, an aligner or a bite splint needs to be worn at night.
- Receding gums or unstable dentition (bite pattern) are present.

- Pain or discomfort in the jaw joint.
- If there is any acute inflammation or there are wounds/injuries present in the mouth or dental area (although in this case it is only necessary to temporarily suspend use, i.e. while the impairment is active).
- If the jaw has not yet finished growing. For this reason this device is not recommended for persons under 18 years of age.
- If you have any loose or unstable bridges.
- If you suffer from chronic asthma or epilepsy.

It is generally recommended that patients have their dentition (bite) and jaw checked by their doctor prior to their first use.

3. Who is mandibular advancement device therapy for?

Therapy using a mandibular advancement device is aimed at two groups of patients. Firstly, patients with position-dependent, mild to moderate obstructive sleep apnoea. For these people, sleeping on their backs causes the tongue muscle to fall backwards into the throat, obstructing the airway. This results in interruptions in breathing.

Secondly, the therapy is intended for patients with simple, position-dependent snoring (also known as primary snoring). With this type of presentation, also known as tongue snoring, the tongue slackens when the patient sleeps on their back, falls backwards into the throat, and narrows the upper airways. Air turbulence at the points of constriction causes the tissue of the throat to vibrate, resulting in a snoring noise.

4. Product features and expected clinical benefits of therapy

Performance features

The mandibular advancement device works mechanically to prevent

narrowing or complete obstruction of the upper airway when sleeping on the back. The mouthpiece shifts the lower jaw forward and thus applies a constant, slight tension to the tongue musculature attached to the lower jaw. This means that the tongue muscle, even when in a state of full relaxation, is prevented from slackening to such an extent that it narrows or closes off the upper airways in the throat area.

Expected clinical benefits of the therapy

The reduction or complete elimination of snoring or interruptions in breathing triggered by sleeping position – meaning no disruptions to breathing during sleep.

5. Risks and unwanted effects

- Tooth sensitivity: For the first minute or so after removing the mouthpiece in the morning, there may be a slight increase in tooth sensitivity. This subsides after a few minutes.
- Pain at the base of the ear: During the first few nights after starting the therapy, painful tension may occur at the base of the ear. This usually disappears after a few days. If the pain persists, please reduce the degree of advancement on the mouthpiece so as to reduce the tension on the jaw. If this does not help, please discontinue the therapy and talk to your doctor.
- Increased salivation: In the first days of use, heavy salivation sometimes occurs when the mouthpiece is inserted in the mouth. The body first needs to get used to having a „dry foreign object“ in the mouth. This can be remedied by wetting the mouthpiece with water before inserting it (= rinse briefly under clear water prior to insertion). After a few days, salivation returns to normal.
- Dry mouth: In rare cases, dry mouth can occur as a result of

mouthpiece therapy. Adhesive tablets for dry mouth (e.g. XyliMelts) can relieve this.

- Nausea: In very rare cases, nausea may develop while using the mouthpiece. In this case, please stop the therapy and talk to your doctor.
- Change in dental occlusion (the contact position of the teeth of the upper and lower jaw when the jaw is closed): It is possible for changes to occur in the contact position of the teeth of the upper and lower jaw as a result of mouthpiece therapy. If you notice a change in the resting position of your lower jaw after wearing the mandibular advancement device (i.e. the lower jaw does not return to its original position), please discontinue the therapy and talk to your dentist.
- Spitting out: If the mouthpiece has not been fitted properly and is too loosely adhered to the dentition (bite), it can be that the mouthpiece is unconsciously spat out at night. As a result, the sleep apnoea or snoring may start up again until the mouthpiece is reinserted.
- Minor burns to the skin if the mouthpiece is not fitted correctly: If the detailed instructions for fitting the mandibular advancement device are not followed, under certain circumstances minor burns can occur to the skin. Fitting requires using hot water at approx. 75° C.
- Swallowing: Accidental swallowing of the mouthpiece is not possible. The mouthpiece is so large that it cannot be easily swallowed.

Should the unwelcome side effects persist, stop wearing the mouthpiece and seek medical advice.

6. Preliminary steps – fitting

Before the first use, the mouthpiece needs to be fitted to the patient's bite.

General information

The water used to fit the mouthpiece must be at a temperature of at least 75° C. If the temperature of the water is too low, the material of the mouthpiece will not soften enough in the water and will not be easily moulded. If the water temperature is too high, the material of the mouthpiece becomes very warm. This can make the fitting uncomfortable for the teeth.

Fitting

1. Take the upper part of the somnipax guard (standard / large) (= with 1 hook, front) and stand in front of a mirror. Now practise inserting the mouthpiece into your mouth in front of the mirror and carefully bite into the groove of the mould. Don't worry if the blank mould feels large and bulky – the mouthpiece will fit in your mouth much better after the fitting process and will barely be noticeable.

2. Place the thermometer (the black disc with the green centre) in a pan filled with about 8 cm of cold water. Heat the water. Stop heating when the temperature indicator in the middle of the thermometer changes from green to red (optimum temperature) **a**.

CAUTION: The water must not reach boiling temperature .

3. When the water has reached the optimum temperature (= the temperature indicator on the disc has changed from green to red), place the upper piece of your somnipax guard (standard / large) into the hot water for 3 minutes. For the following steps, please note that the fitting must be done within 10–15 seconds of removing the mouthpiece from

the hot water. Otherwise, there is a risk that the material will cool down too much and it will no longer be possible to shape it.

4. Now remove the upper jaw part from the water using a spoon and shake off the hot water. Ensure all hot water is removed.

CAUTION: Before inserting, always perform a check (e.g. by holding it briefly against the wrist) to ensure the mouthpiece will not cause burns in the mouth.

Now, standing in front of the mirror, insert the upper part of the mouthpiece onto the upper jaw area using a swift movement. Please align the mouthpiece so that the white hook on the upper jaw part is positioned exactly in the centre between the two front incisors . Now press the upper jaw part upwards into your upper jaw area on both sides using both thumbs **b**.

5. Hold the mouthpiece in your mouth with one hand and insert the bite plate swiftly with the other hand. Now bite down hard to get a good imprint of your teeth **c**. After that, press the inserted mouthpiece against the side of your teeth with your fingers. Reduce the bite pressure slightly and let the mouthpiece cool and harden in your mouth for 2 minutes. Keep the jaw as still as possible during this time. The fitting of the upper jaw part is now complete.

6. Please keep the upper jaw part inserted in your mouth during the fitting of the lower jaw part, or reinsert it. The lower jaw part (= the part of the mouthpiece with 2 hooks on the side at the back) is fitted in exactly the same way as the upper jaw piece. I.e. heat the water to the optimal temperature (= until the temperature indicator on the disc has changed from green to red), then leave the lower jaw piece in the hot water for 3 min.

Remove the lower jaw part from the hot water using a spoon and shake

it to completely remove any hot water **d**. Check once again that the MAD will not burn your mouth before inserting it. Insert the lower jaw part into the mouth so it aligns with the upper jaw part, i.e. their two vertical markings must be level. Press down firmly with both index fingers and press the two mouthpiece parts together by biting down with increasing pressure **e**. Then reduce the pressure slightly and allow the two mouthpiece parts to cool and harden in the mouth for at least 2 minutes. After that, rinse the lower part in clear, cold water. The lower jaw part is now ready.

7. The two mouthpiece parts now need to be connected using the tension bands. Using scissors or a knife, cut band number “4” (or possibly “5” depending on the model, see below) from the band holder **f** and insert it. First guide the hook on the upper jaw part (= the mouthpiece parts with only one hook at the front) into the central loop of the band. To do this, turn the mouthpiece 90 degrees so that the long holding pin can simply slide over the (elongated) band loop. Then turn the mouthpiece back by 90 degrees again – the band is now firmly fixed and can no longer slide off the hook. Repeat the process on the lower jaw part using the two elongated loops on the ends of the band (turn 90 degrees > let it slide over the elongated loop > then turn back 90 degrees > done). The mouthpiece is now ready for use!

NOTES ON THE CORRECT BAND LENGTH / CORRECT DEGREE OF ADVANCEMENT:

There are large number markings on the bands of the somnipax guard (standard / large) mouth guards. The bands cause the following advancement of the lower jaw:

somnipax guard standard	
Band no.	Advancement (mm)
Number 1	0,0
Number 2	3,0
Number 3	4,5
Number 4	6,0
Number 5	7,5
Number 6	9,0
Number 7	10,5
Number 8	12,0

somnipax guard large	
Band no.	Advancement (mm)
Number 1	-1
Number 2	0,5
Number 3	2,0
Number 4	3,5
Number 5	5,0
Number 6	6,5
Number 7	8,0
Number 8	9,5

- **somnipax guard standard:** you should use band length 4 for the first few nights after adjustment (this has an advancement of 6 mm).
- **somnipax guard large:** you should use band length 5 for the first few nights after adjustment (this has an advancement of 5 mm).
- Shorter/longer bands can then be used after the first few nights as required. The shorter the band, the larger the advancement.
 - If band 4 is too tight, for example (= there's an unpleasant pulling feeling in the jaw), a longer band should be used, i.e. a band with a lower number (1, 2 or 3).
 - If you still experience breathing stops or snoring when using band 4, a shorter band should be used, i.e. to increase the advancement (bands 5, 6, 7 or 8).
- An advancement of 4-7 mm is usually adequate for most users (i.e. the corresponding bands are the most commonly used). If the mouth

guard is being used to treat obstructive sleep apnoea, the handling physician should set the correct advancement or check the setting you have chosen.

General information about fitting:

- The better the quality of the dental impression, the better the hold of the mouthpiece.
- The fitting can be repeated several times if necessary, as the material is thermoplastic and can be heated multiple times while retaining its properties.
- Any loose edges created during the fitting can be carefully removed with a sharp knife.

7. How do I use the mandibular advancement device correctly?

Before inserting the mouthpiece, i.e. before going to bed, clean your nasal passages, for example using a neti pot. Give the clean mouthpiece a quick rinse under cold water and insert it into your mouth. Experience tells us that the easiest way to do this is if you first bite into the upper mouthpiece with your upper jaw teeth and then pull your lower jaw forward and bite into the lower mouthpiece with your lower jaw teeth.

Please bear in mind that it may take a few nights to fully get used to the mouthpiece and for the sensation of a foreign object in your mouth while wearing it to subside.

The mouthpiece should be rinsed with cold water each time you remove it from your mouth. Please follow the cleaning instructions.

8. How do I determine if the mouthpiece is effective or being used correctly?

Habitual snoring

The easiest way is to ask your sleeping partner if you're still snoring

when sleeping on your back (please note that even with effective anti-snoring therapy, a certain residual noise level may persist due to natural breathing sounds).

If you sleep by yourself, you can use a smartphone app to record your snoring. This can help you check whether the therapy is working effectively. A dry mouth in the morning is also an indication – albeit a less decisive one – that breathing was taking place through an open mouth and snoring probably continued.

Obstructive sleep apnoea

Here too, you can ask your sleeping partner whether interruptions in breathing are still occurring when you sleep on your back. However, as treatment for obstructive sleep apnoea is always carried out under the supervision of a doctor, this assessment alone will not be enough. Your doctor will also need to determine whether the therapy is effective (sleep endoscopy examination / polygraphy).

9. How do I deal with impaired functioning of the mouthpiece? Nighttime spitting out / loss

If the mouthpiece has not been fitted properly, i.e. the material of the mouthpiece does not adhere closely enough to the surface of the teeth, it may be spat out during use. As a result, the advancement effect is lost and the therapy is not effective. In this event, clean the mouthpiece and reinsert it. If it still does not form a tight fit, the fitting process must be repeated.

Points of pressure from the mouthpiece

It may be that after a successful fitting the mouthpiece feels too large or presses onto the gums in some places (this issue may occasionally occur on the surface of the teeth). You feel less willing to wear the

mouthpiece as a result. Check the fit of the mouthpiece, and at the pressure points on the gums carefully trim the edges with a sharp knife or cutter. Check that it fits correctly.

Deformation due to heat

If the mouthpiece is exposed to excessive heat over a long period of time or is stored in too hot an environment (especially in direct sunlight), it can become deformed. It then becomes unusable and must be replaced.

10. When should I seek medical advice?

If any of the following apply, we recommend talking to your doctor before using the mouthpiece.

- Symptoms relating to the temporomandibular joint (TMJ) (► dentist / orthodontist)
- Fewer than 10 firmly seated teeth per jaw (upper / lower jaw) or unstable dentition (► dentist)
- Patient uses a full or partial removable denture or nighttime braces in the mouth or jaw area (► dentist)
- Receding gums, acute inflammation or wounds/injuries in the mouth or tooth area (► dentist)
- Protruding lower jaw/pronounced overbite (► dentist/orthodontist)
- Technical issues with fitting the mouthpiece (► here we recommend asking the dentist or ENT doctor responsible to take over or assist with the fitting – if necessary, with a contribution to the costs)

11. Correct cleaning, care and storage

Cleaning and care

For everyday use, rinse your mouthpiece with cold water only and store it in a cool, dry place in the plastic box provided. At regular intervals,

and at least once a week, the mouthpiece should be cleaned with a cleaning agent formulated for mouthpieces.

To do this, dilute the cleaning agent with lukewarm water as specified by the manufacturer and clean the mouthpiece briskly with a soft toothbrush. Then rinse the mouthpiece thoroughly with cold water.

Please do not soak the mouthpiece in the cleaning solution for an extended period of time. The surface of the mouthpiece is delicate and could otherwise be damaged by the solution.

Alternatively, the mouthpiece can be cleaned using cleaning tablets („tabs“) for mandibular advancement devices. Please do not use off-the-shelf commercial cleaning tablets designed for dental or maxillary prostheses or cleaners for braces. These are too aggressive for the soft surface of the mouthpiece and will corrode or damage it.

Storage

Please store the mouthpiece in a cool, dry place in the storage box provided, out of reach of children and pets. The lifetime of the mouthpiece is usually 8–18 months. Extra stress (severe bruxism, high mouth acidity...) can significantly shorten its lifetime.

12. How can I determine if the mouthpiece is still ok to use?

If the surface of the mouthpiece is cracked, broken, or otherwise damaged, it needs to be replaced. Please dispose of mouthpieces that are damaged or need to be replaced in your household waste.

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