**Gebrauchsinformation/Manual/
Notice d'utilisation/Indicazioni d'uso/
Instrucciones para el uso/Handleiding/
Käyttöohje/Instruktioner för användning**

DEUTSCH	3
ENGLISH	9
FRANÇAIS	15
ITALIANO	21
ESPAÑOL	28
NEDERLANDS	34
SUOMI	40
SVENSKA	46

Lieferumfang/Delivery includes/Inclus/Ambito di consegna/ Alcance de la entrega/Omvang van de levering/ Pakkauksen sisältö/Omfattning av leveransen

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Manual ENGLISH

1. What is somnipax belt and what is it used for?

somnipax belt is an active positional therapy device used for the relief of simple snoring caused by sleeping on the back (also known as tongue snoring), as well as for mild or moderate positional sleep apnoea. somnipax belt is a class 1 medical device, Regulation (EU) 2017/745.

2. Intended use and scope of application

The intended use is positional therapy for snoring related to sleeping on the back, as well as for mild and moderate sleep apnoea. If you are in a back-sleeping position that promotes snoring / sleep apnoea, the device vibrates until you change position onto your side or front.

It should not be used:

- by wearers of cardiac pacemakers,
- during convalescence following intervertebral disc surgery,
- if you suffer from severe spinal complaints,
- positional vertigo,
- position-independent central sleep apnoea syndrome,
- by persons under 18 years of age.

 Should you have any doubts (e.g. because of possible interactions between medication and the positional therapy), please talk to your doctor.

3. For whom is the therapy with somnipax belt intended?

Therapy with somnipax belt is aimed in the first instance at patients who suffer from positional, mild or moderate sleep apnoea. In this group of

patients, sleep apnoea is caused by the slackened tongue base muscle slipping into the throat during sleep and thereby temporarily closing the airways. This leads to interruptions in breathing. If you are in a back-sleeping position that promotes snoring / sleep apnoea, the device vibrates until you change position onto your side or front. This prevents the tongue base muscle from slipping into the throat.

In non-pathological tongue snoring, the slackened tongue base muscle again slips into the throat during sleep, but only narrows the throat rather than obstructing it entirely. The narrowing of the throat causes the snoring sounds. somnipax belt prevents this narrowing of the throat by vibrating to signal a back-sleeping position (which promotes snoring / sleep apnoea) until you change position onto your side or front.

Contraindications:

Wearers of cardiac pacemakers must not use this device.

4. Product features and expected clinical benefits of therapy

Product features

somnipax belt vibrates to signal that you are in a back-sleeping position (which promotes snoring / sleep apnoea) and does so until you change position onto your side or front. The slackened tongue base muscle can then no longer obstruct or significantly narrow the throat.

Expected clinical benefits of the therapy

Reduction or complete prevention of interruptions in breathing caused by obstruction of the throat during sleep (positional, mild and moderate sleep apnoea). Reduction or complete prevention of tongue snoring.

5. Risks and unwanted effects

The somnipax belt requires an adjustment period of approximately 1-2

weeks. During the initial nights in particular, sleep may be disturbed by the vibration. The belt should be worn snugly against the body to ensure effective transmission of the vibration. If this is not observed, more vibrations may need to be triggered, resulting in increased noise which could disturb your partner's sleep. Should you experience any discomfort in your spine or musculoskeletal system (particularly in the shoulder, neck or throat area) during use, please consult your GP before continuing use.

Warning and safety instructions

Please read through the instruction manual. Failure to follow the terms of use and instructions contained in this instruction manual may affect the safety and effectiveness of the device. The manufacturer and distributor do not bear liability if the device has been modified or repaired by unauthorised persons or if the device has been subject to any other form of technical intervention by unauthorised persons.

- Please also read the instruction manual before using the product. The instruction manual is an integral part of the device. Please store it in the same place as the device. Only use the device for its intended purpose.
- If you experience daytime sleepiness, you may be suffering from sleep apnoea; in this case, please be sure to visit your doctor for an in-depth assessment.
- When opening the packaging, check that the device is undamaged, paying particular attention to the plastic parts.
- Only use alkaline batteries of the type specified in this instruction manual. Please remove the batteries from the device if it will be left unused for an extended period.

- Wear the device over undergarments or nightwear.
- Do not use the device in the presence of flammable anaesthetic mixtures containing air, oxygen, or nitrous oxide.
- This medical device requires special safety precautions with regard to electromagnetic compatibility (EMC) and must not be used near portable or mobile RF communication devices (mobile phones, radios, etc.) as they may interfere with its operation.
- Do not touch the device with wet hands, and protect it from contact with liquids.
- In the event of breakage or malfunction, stop using the device immediately and contact your retailer.
- Children and/or people with learning disabilities may only use the device under the strict supervision of a responsible adult.
- Repairs may only be carried out by the manufacturer's repair service or a service centre authorised by the manufacturer with the use of original spare parts. Making any modifications to the device without prior authorisation from the manufacturer is prohibited. With the exception of replacing the battery, no electrical and/or mechanical components of the device are designed to be repaired or modified by the user.
- Please take note of the temperature ranges indicated in the technical data.
- If you show any signs of hypersensitivity or allergic reactions to the device materials, stop using the device and consult a medical specialist. Please also notify the manufacturer or retailer.

6. Preliminary steps

Please refer here to the information in 7.

7. How do I use somnipax belt correctly?

Fitting the belt:

- Remove the battery protection film by pulling the red tab.
- Insert the device into the belt pouch.
- Place the belt around your chest with gentle tension.

How it works and how to use it

To switch on the device:

- Press the on/off button briefly (less than one second). Two short vibrations tell you that the system is active. You can relax and lie down on your side or front.
- If you roll over into a back-sleeping position, after 10 to 15 seconds the system will begin to produce vibrations at regular intervals until you move back into a side- or front-sleeping position.

To switch off the device:

- Press the on/off button for 2 to 3 seconds. A long vibration indicates that the device has switched off.

Adjust the intensity of vibration:

- With the device switched off, hold down the on/off button for at least 5 seconds. The device will produce vibrations of varying intensity. Release the button once the vibrations have reached a comfortable level. The device has now saved your vibration setting and remains switched off and ready to start.

8. How do I determine whether somnipax belt is working and/or being used correctly?

Habitual snoring

The easiest way to do this is to ask your sleeping partner whether you are now sleeping on your side or stomach and whether you are continuing to

snore while you sleep (please note that even with effective anti-snoring therapy, a certain residual noise level can persist due to natural breathing noises).

If you sleep alone, smartphone apps can be used to record snoring. These can help you check whether the therapy is working effectively. A dry mouth in the morning is also an indicator – albeit a weaker one – that breathing took place through an open mouth and snoring probably continued.

Obstructive sleep apnoea

Here too, your sleeping partner can be asked whether you are now sleeping on your side or on your stomach and whether interruptions in breathing are still occurring. Since treatment for obstructive sleep apnoea is always under the supervision of a doctor, however, this assessment will not be enough. Your doctor will also need to determine whether the therapy is effective (sleep endoscopy examination / polygraphy).

9. What do I do if the function of the positional therapy belt is impaired?

The battery is largely or completely discharged (red indicator light is on):

1. Remove the device from the belt. Unscrew* the two screws on the back of the housing. Remove the bottom of the housing.
2. Take out the empty batteries and press the On/Off button several times. Insert two new AAA 1.5V batteries.
3. When inserting the batteries, the system will emit a vibration and the GREEN INDICATOR light (new batteries) will illuminate for approximately 1 second**. The vibration is set to maximum strength by default.

4. Close the device in the correct orientation without using force; if necessary, rotate it half a turn and tighten both screws*. Insert the housing back into the belt pouch.

**Note: Use only a PH1 (Phillips head) screwdriver.*

***Note: If the device does not respond to the insertion of new batteries, check that the batteries are inserted with correct polarity. If the problem persists, please verify that the inserted batteries are new and charged.*

The fabric belt has significant nicks or tears or the device itself is cracked:

If the fabric belt displays significant nicks or even tears or the device itself is cracked, please stop using somnipax belt immediately and replace the damaged part. Slight imperfections on the fabric belt (e.g. loosened seam or loose threads) will not generally affect its function.

10. When should I seek medical advice?

If you show any signs of hypersensitivity or allergic reactions to the device materials, stop using the device and consult a medical specialist. Please also notify the manufacturer or retailer.

11. Correct cleaning, care and storage

Cleaning and care

Clean the device with a damp cloth. Do not rinse under running water. The fabric belt can be washed by hand as per the instructions on the wash label.

Protection against water: Not protected!

Storage

Please store your somnipax belt in a dry place away from direct sunlight.

Permissible storage temperatures for somnipax belt during storage and transport are between -20 and +50 degrees Celsius (at a humidity of 0% to 80%). During use, permissible temperatures are between 10 and 30 degrees Celsius (at a humidity of 0% to 80%).

12. How do I determine if my somnipax belt is still ok to use?

Please always check that the device is still working correctly before use: When the system is switched on, the green indicator light will flash at regular intervals (approx. every 5 seconds). Then also briefly pressing the on/off button causes a short vibration. Note: When the system is switched off, the green indicator light flashes every 30 seconds.

Battery display:

If the red indicator light is on instead of the green indicator light, this means the batteries will soon need to be changed. If the batteries are not replaced in time, the device will not switch on and three consecutive pulses of vibration and flashes of the red indicator light will be shown.

Automatic switch-off:

To prevent the batteries running down prematurely, the device switches itself off after approximately 14 hours.

If the device is cracked or otherwise damaged, it must be replaced. The same applies if the fabric belt is damaged or defective. The fabric belt can be disposed of with your household waste. The device itself must be disposed of in accordance with the regulations in effect for electrical appliances.

13. Technical data

Power supply	2 alkaline batteries AAA 1.5 V
--------------	--------------------------------

Battery life	2 to 6 months*
Classification	Class 1 medical device in accordance with (EU) 2017/745.
Vibration function	Adjustable intensity of vibration, 7 levels
Display elements	One red indicator light and one green indicator light
Belt size	circumference 70 cm to 125 cm
Weight	26 g (without batteries)
Dimensions electronic positional therapy device	75.5 mm x 37.5 mm x 14.5 mm**
Operating conditions	room temperature: 10÷37 °C humidity: 0÷80% rF (free of condensation)
Storage and transport conditions	room temperature: -20÷50 °C humidity: 0÷80% rF (free of condensation)
Applied part type	BF

**Calculated using brand new, fully charged, high quality batteries available as standard on the market. The minimum battery life of two months is guaranteed with a total maximum vibration time of one hour per night.*

*** Device without fabric belt.*

14. Information on environmentally friendly disposal of electrical appliances

1. Disposal of waste electrical and electronic equipment (WEEE)

The symbol of the 'crossed-out wheelie bin' means that you are legally required to dispose of these appliances separately from your residential waste. Disposing of them via household waste, such as in the recycling bin or general waste bin, is prohibited. Avoid waste being misdirected by disposing of it correctly at special collection and return points.

2. Removal of batteries and bulbs

This product contains 2 AAA 1.5 V alkaline batteries which can be removed from the end-of-life device without causing any damage. These must be removed and disposed of as batteries separately before disposing of the device.

3. Options for returning end-of-life devices

Owners of end-of-life devices can return them free of charge within the range of options provided by public waste management authorities for the return or collection of WEEE, in order to ensure proper disposal of the end-of-life device.

4. Data protection

We would like to point out to all end users of WEEE that you are responsible for deleting personal data on the end-of-life device to be disposed of.

5. WEEE registration number

We are registered as a manufacturer of electrical and/or electronic equipment under the registration number DE10291720 with the Stiftung Elektro-Altgeräte (stiftung ear) register, Nordostpark 72, 90411 Nuremberg, Germany.

6. Collection and recycling quotas

In accordance with the WEEE Directive, EU member states are obliged to collect data on WEEE and to transmit this data to the European Commission. Further information is available on the website of the Federal Ministry for the Environment, Nature Conservation, Nuclear Safety and Consumer Protection.

15. Information on environmentally friendly disposal of batteries

As the battery end user, you are legally obliged to return used batteries. Batteries and rechargeable batteries must not be disposed of with household waste. After use, you can return batteries free of charge at a retail outlet.

Meaning of the symbols on the batteries:

- The crossed-out wheelie bin means that the product must not be disposed of with household waste.
- Hg = battery contains more than 0.0005% mercury by weight; Cd = battery contains more than 0.002% cadmium by weight; Pb = battery contains more than 0.004% lead by weight.

Avoid waste by using batteries with the longest possible battery life or rechargeable batteries („accumulators“). Please avoid polluting our environment with waste and do not throw batteries away carelessly or leave them lying around. Before disposing of a battery according to regulations, please check the options available for preparing used batteries for reuse. Please be aware of the possible harmful effects of the substances contained in batteries on the environment and human health, particularly the risks associated with handling batteries containing lithium. Separate collection and recycling of used batteries is intended to help minimise the

impact of batteries on the environment and human health.

16. Required information on electrical safety

Device class	I
Applicable standard	IEC 60601-1: 2005 + A1:2012
Applied part type	BF
Sterilisation	None

17. Required information on electromagnetic compatibility

Classification: The device has been designed to comply with the requirements of the IEC 60601-1-2:2015 standard on electromagnetic compatibility.

However, certain materials may inadvertently emit strong RF signals which may interfere with the device

It is recommended that the device be kept away from other electrical equipment (in accordance with **section 7.9.2.2 IEC 60601-1-2**).

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Table of tests

Test Designation	Severity	Verdict	Basic Standards/Comments
Harmonics current emission	-	N/A	EUT powered through battery
Measurement of voltage fluctuation and flicker	-	N/A	EUT powered through battery
Measurement of discontinuous disturbances	-	N/A	EUT powered through battery
Measurement of radiated disturbances			EN 55011
-Front side - Patient on right position mode	B	PASS	
-Lateral side - Patient on right position mode	B	PASS	
-Lateral side 1 - patient alarm mode	B	PASS	
-Lateral side 2 - patient alarm modew	B	PASS	

Electrostatic discharges immunity			EN 61000-4-2
-INDIRECT 8kV (VERT.)	I.3 (60601-1-2)	PASS	
-INDIRECT 8kV (HORIZ.)	I.3 (60601-1-2)	PASS	
-DIRECT 8kV	I.3 (60601-1-2)	N/A	
-AIR 2kV	I.3 (60601-1-2)	PASS	
-AIR 4kV	I.3 (60601-1-2)	PASS	
-AIR 8kV	I.3 (60601-1-2)	PASS	
-AIR 15kV	I.3 (60601-1-2)	PASS	

Radiated, radio-frequency, electromagnetic field immunity -Front side - Patient on right position mode -Front side - Patient on right position mode (spot frequencies) -Front side - Patient on right position mode (spot frequencies)>1GHz -Front side - Patient on right position mode (spot frequencies)>2.7GHz -Lateral side - Patient on right position mode -Lateral side - Patient on right position mode (spot frequencies) -Lateral side - Patient on right position mode (spot frequencies)>1GHz -Lateral side - Patient on right position mode (spot frequencies)>2.7GHz -Lateral side 1-patient alarm mode -Lateral side 1-patient alarm mode (spot frequencies) -Lateral side 1-patient alarm mode (spot frequencies)>1GHz -Lateral side 1-patient alarm mode (spot frequencies)>2.7GHz -Lateral side 2-patient alarm mode -Lateral side 2-patient alarm mode (spot frequencies) -Lateral side 2-patient alarm mode (spot frequencies)>1GHz -Lateral side 2-patient alarm mode (spot frequencies)>2.7GHz	I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2) I.3 (60601-1-2)	PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS PASS	EN 61000-4-3
--	---	---	---------------------

Electrical fast transient/burst immunity	I.3 (60601-1-2)	N/A	EUT powered through battery and without cables
Surge immunity	I.3 (60601-1-2)	N/A	EUT powered through battery
Conducted disturbances induced by radio-frequency fields immunity	I.3 (60601-1-2)	N/A	EUT powered through battery and without cables
Power frequency magnetic field immunity	I.3 (60601-1-2)	N/A	EUT has no sensitive components
Voltage dips and short interruptions immunity	I.3 (60601-1-2)	N/A	EUT powered through battery
Measuring electromagnetic fields for comparison with reference levels -1Hz - 400kHz -100kHz - 3GHz	6.25µT 27.5V/m	PASS PASS	EN 62311:2008

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EN: Should any serious incident occur in connection with the use of the somnipax belt, please report it to the relevant authority in your country of residence and to the manufacturer.

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